

## MEDOR: special application

**Objective:** H<sub>2</sub>S and THT analysis in a matrix with 10% of H<sub>2</sub>.

The analyser is calibrated on DMS compound (internal standard) in middle amplification from DMS at 33.7 ppb (+/-8%)

The response factors are (for the same amplification):

R<sub>f</sub> (DMS) = 1

R<sub>f</sub> (H<sub>2</sub>S) = 2.22

R<sub>f</sub> (THT) = 0.35

The **H<sub>2</sub>S** compound comes from a permeation tube heated and swept by zero air.

The **THT** compound comes from a cylinder diluted in zero air.

### Analytical conditions

Packed columns at 40°C – 1/8" Teflon – length: 190 cm + 40 cm

Carrier gas: N<sub>2</sub> at ≈ 5 ml/min – Pressure: 652 hPa


Loop: 1 ml - Sampling flow: - 100 ml/min

Analysis cycle: 600 s – Acquisition: 450 s – Amplification: high (3)

Detector: wet cell.

### 1. H<sub>2</sub>S at 26,8 ppb and THT à 19.6 ppb (+/- 10%) in a matrix with 10.2% H<sub>2</sub>

Chromatogram :



Information and operating conditions

Analyser :

Serial Number : #44520111

Owner : AIR LIQUIDE

Location : JOUY EN JOSAS

Operating conditions :

Description : Sample analysis-A=3

Method Name : SAMPLE3

Substances Table Name : SOUFRE3

Sampling :

Tube Number : 1

Duration : 0 s

Volume : 1,0 ml

Date : 26/01/2011 03:29:00

Detector :

Amplification : 3-High

Sample Rate : 6 per second


Peak List

| Substance | Result | Unit | Start  | R.Time | Max      | Stop   | Area     | Type |
|-----------|--------|------|--------|--------|----------|--------|----------|------|
| H2S       | 26,76  | PPB  | 12,17  | 22,83  | 20239,00 | 50,00  | 23233,00 | ST_E |
| THT       | 18,97  | PPB  | 227,50 | 261,83 | 25934,00 | 325,17 | 2,7034E5 | ST_E |

## 2. 24 hours trend on standard gas



Trend on retention times in seconds



Trend on concentrations in ppb