

Chroma CO

A gas analyser using flame detection with a catalytic methaniser module to analyse CO / CH₄ and CO₂ for air or gas monitoring



Integrated computer version (5U high) : C11022

Industrial hygiene / Process / Environment / Emissions

ambient air:

CO / CH₄ and CO₂ in ambient air
Toxic compounds monitoring: **HCHO** and **Acetaldehyde** in option

In air (CO = around 300 ppb) (CO₂ = 380 ppm and CH₄ = 1800 ppb)

Process / Pure gas:

Methane / CO / CO₂ in pure gas
Monitoring of CO in pure N₂O
Quality control of manufacturing process

* in option

Chromatotec® is specialized in VOC, Sulfur and permanent gases analysis at trace and ultra trace levels (ppm, ppb, ppt).
Please visit our website for more details.

Updated: April 2021

Chroma CO



A gas analyser using flame detection with a catalytic methaniser module to analyse CO / CH₄ and CO₂ for air or gas monitoring

Principle:

A 19» inch rack mountable chassis, 5U high. It has the ability to be used as a fully integrated system for automatic sample analysis as well as a transportable unit for use on various sites. Our application flexibility allows the chromaCO start up to meet the customer's specifications.

The CHROMA CO is an **isothermal industrial gas chromatograph** (isothermal column) dedicated to the analysis of CO/CO₂/CH₄/ in air or industrial hygiene or pure gas like N₂. As a fully automatic instrument, it allows **unattended operation for commissioning 24/24**. The injection is made via a heated valve. CO CH₄ CO₂ are separated by the column and CO and CO₂ are hydrogenated in CH₄ before the FID.

- A valve regulates pure carrier gas.
- The instrument has a flame ionization detector FID.
- The detection limit is in ppb in standard.
- It is equipped with an internal microprocessor which runs all the analytical events (temperature, pressure injection, and detection and concentration calculations).

The Vistachrom software enables the user to visualize and store data on a computer. Furthermore it provides comfortable utilities to recalculate, calibrate and export data and to set-up measurement. The airmo TREND software allows the calculation of retention time, area, mass or concentration profiles.

CALCULATION AVERAGE can be made for better precision

Options:

On-line results are transmitted via :

- NMTHC (+ THC by calcul = NMTHC + CH₄)
- HCHO application
- A Modbus / Jbus or BH protocol.
- Analog output 4-20 mA or 0-10 V.
- **Automatic validation**
- **Generator of purified air (air with THC in ppt).**
- **Oven gradient temperature**
- 24 V power supply.
- Multiple stream selector (2 to 6)
- Pressurised housing like Exp for installation in area 1..
- Available in version integrated computer – 5U (including ref. XXX022 Electronic and LCD display).
- **Gas generator H₂/air/N₂**

Product technical specifications:

CO/CH₄/HCHO and CO₂ analysis:

- For information in air : CO = 300 ppb : CO₂ = 380 ppm and CH₄ = 1800 ppb

Detection limit for CO CH₄ CO₂:

ambient air:

- 50 ppb with H₂ (5.5)

pure gas:

- 50 ppb with N₂ (5.5)
- 5 ppb with N₂ purifier

Detection range:

ambient air:

- 0 - 10 ppm (CO/CH₄) & 0 - 1000 ppm (CO₂)

pure gas:

- 0 - 10/20/50ppm

Relative Standard Deviation at middle range :

- < 0.3 % over 48h (RT)
- < 2 % over 48 h (Conc.)

Results:

- Data storage
- 4-20mA current output (option)
- Module for 4 output 4-20 mA current output (option)
- Modbus / Jbus or MGS1 communication protocol (option)

Cycle time:

- 5 minutes for CH₄/CO₂
- 10 minutes for CH₄ CO CO

Gas supply:

- **N₂ as carrier gas for better sensitivity**
- H₂ (FID and carrier gas): 30 ml/min (inlet 2 bars ; 1/16» swagelock)
- Air (FID): 180 ml/min (inlet 3 bars; 1/8» swagelock)
- Sample inlet (vacuum pump) ¼» swagelock
- Pneumatic valve 90ml/commutation

Sample loop:

- Injection loop 0.5 or 1 or 2 ml ,depending of the range.

Power supply:

- Main: (230V / 115V 50 Hz/60Hz)
- Battery: 24V (option)

Electrical consumption:

- 250 VA

Dimensions and weight:

- Rack 482 mm (19"), Height 177mm=4U, Depth 600mm 5U for C11022
- Weight: 20kg

To order:

chroma CO (4U)
Needs a supervisor
Supervisor (5U)

Model:

C11000
XXX022

Chromatotec® is continuously improving its products, therefore these specifications are subject to change without notice

To contact us: info@chromatotec.com



NORTH AMERICA

CHROMATOTEC Inc.
18333 Egret Bay Blvd, Suite 270,
Houston TX 77058 - USA
Phone: +1 (281) 335 4944
Fax: +1 (281) 335 4943

EUROPE

AIRMOTEC AG SAS
15 rue d'Artiguelongue
33240 Val de Virvée, FRANCE
Phone: +33 (0) 557 940 626
Fax: +33 (0) 557 940 620

ASIA

Chromatotec Trading (Beijing) Co., Ltd.
Room 1806, Building 1,
Wanda Plaza, No.93, Jianguo Avenue,
Chaoyang District,
Beijing 100022 - CHINA
Phone : +86 (0) 105 960 3283