

Vistachrom 1.4.9

Changes and functionalities

Chromatotec®

Contents

- Vistachrom updates

What are the changes between 1.4.7 and 1.4.9 ?

- The RealTime DataBase

How does it work ? What is the interest ?

- The Equipment Manager

What is it ? What is it used for ?

- Vistachrom possibilities to meet your requirements

1. Vistachrom updates

(1.4.7 → 1.4.9)

What is Vistachrom ?

- A software suite, but also the main application program of that suite
 - Main application
Vistachrom
 - Linked modules
EquipmentManager, MathModule, MJBus driver, MGS1 driver
 - Standalone applications
MethodManager, PeakViewer, ServiceGC, UnitManager

What is Vistachrom role ?

- Monitoring one or more analyzers
 - Each analyzer is in real time controlled with a CPU board.
Detector signal processing, heating and valves monitoring, ...
 - Vistachrom communicates with those boards: monitors them, sends and retrieves data.
 - Vistachrom is in charge of analysis configuration and data post-processing.
Configuration updating, error checking, data saving, output configuration, ...

What is new with the 1.4.9 version ?

- Internal data management is now centrally-managed : the RealTime DataBase has been created
 - Only one data repository with an user-friendly interface
(UserInterface no longer exists)
 - Modular design : each application program of the suite is independent
- AnalogOutput and DigitalOutput drivers are removed. They are replaced with a more generalist driver offering more opportunities : the EquipmentManager.
- Vistachrom is no longer limited to Windows XP

2. The RealTime DataBase

(RTDB - )


Data in the software suite 1/2

- The real time control of an analyzer is done by the CPU board, the Vistachrom software suite does the supervision.
 - The RealTime DataBase contains analyzer status, detected errors, analysis results
 - No raw acquisition data in it but it can be added. And more widely, all data processed in an application of the suite can be added.
- Data arrangement
 - Before data were “mixed” between applications (1.4.7)
 - Now data are in one reference table : the RealTime DataBase

Data in the software suite 2/2


1.4.7

- Using a shared memory space, applications are interconnected and warn when data is updated.



1.4.9

- An application, the RTDB, centralizes all data. All applications connect to it, subscribe to the data of interest. They are therefore notified by RTDB for these data.



RealTime DataBase specifications

- Each application of the suite connected will subscribe to data interesting it and also update others : it is independent and only works with the RealTime DataBase.
 - Adding a new application to the software suite will be easier
 - Updating applications suite will be less impacting
- The application is launched at Vistachrom startup and it works up to that Vistachrom stops.
- An interface allows you to see all data it contains and the values updating in real time.

The RealTime DataBase

RealTimeDataBase	
Add Delete	
Name	Value
#43590514	
Command_Reg	0
Error_Reg	0
LifeSignal	11896
Parameters	
Ambiante	
Columns	
Pressure	253.487045288086
Temerature	100.2265625
CycleTime	69
LifeTime	11980
PermeationTube	
Results	
TS	
CALIB1	
CALIB2	
COS	
H25-TOS	
BaseSensitivity	5605.73828125
CycleTime	240
SamplingDate	2014-07-03T12:51:00
Substances	
H25	
PeakArea	31263.900390625
Result	4.07130098342896 mg/m³
RetTime	44.7000007629395
Status	0
TOS	
PeakArea	531824.3125
Result	94.8714141845703 mg/m³
RetTime	60
Status	0

- All data of the software suite needing to be exchanged between applications is in this data base.
- Data saved and displayed in a tree structure.
- Displays the added data and updated values in real time.

3. The Equipment Manager

(EM - )


What is EquipmentManager role ?

- Application to communicate with external instruments using the Modbus communication protocol
 - Modbus, Master node (*MJBus driver is a slave node*)
 - RTU, ASCII and TCP compliant
- Replaces AnalogOutput and DigitalOutput drivers because it can control our output modules.
- More widely, it allows to communicate with a lot of industrial devices. And so, it is useful to bring back information and data in the RealTime DataBase.
 - For example, it can communicate with the H2 generator and bring back its pressure in the RTDB

EquipmentManager specifications

- Application configuration “open”
 - First level for customers, distributors, after-sales service, ...
 - Second level for more opportunities but dedicated to the software department
- Two interfaces but little useful for a daily use
 - The first to check communication log and control the links “RTDB data ↔ instrument output”
 - The second to force output or send Modbus requests is dedicated to the maintenance

The EquipmentManager



- Currently used for the 4-20mA and relay outputs
- An example of improved outputs, the 0-4/20mA : thanks to the manager configuration and to the post-processing possibilities of Vistachrom


4. Vistachrom possibilities to meet your requirements

Data outputs proposed

- The device containing Vistachrom can output its data within :
 - 4-20 mA output
 - Relay output
 - Modbus/JBus output (only serial port, RTU-ASCII)
 - MGS1 output (serial port)
- For particular cases, there is the possibility of adding physical gateways and therefore increase the number of proposed choices :
 - Profibus
 - Modbus output for Ethernet (TCP) : *under study*

Data post processing possibilities

- The MathModule allows to retrieve data from the RealTime DataBase, to process them and also to add new ones.
- The MathModule is an opened access to data which offers post processing possibilities on-demand.
 - 0-4/20mA : *Set 0 mA to the output if an error occurred during the analysis*
 - TotalSulfur : *Control analysis and sum defined substances to obtain the total sulfur concentration*
 - Alarm&Stats : *Alarm positioning and statistics calculations*



Conclusion

- New Vistachrom version
 - Modularity and better data management thanks to RealTime DataBase
 - Ease of integration for external devices thanks to EquipmentManager
- Vistachrom offers many opportunities to meet the needs
- But software department is often essential to achieve it



Thank you for your attention !