

AirmoVOC C6-C12 Linearity / Repeatability Test 10, 5, 1 et 0.1 ppb


We have tested our analyser airmoVOC C6-C12 with samples with different concentrations in order to check if the signal varies linearly with the concentration.

We did 3 series of 10 successive analysis with concentration of 10, 5, 1 and 0.1 ppb.

Samples used comes from a standard gas cylinder of 58 compounds. The cylinder is certified at $\pm 2\%$ for each compound and the concentration of the compounds into the bottle are between 98 ppb and 100ppb (see appendix) so we need to dilute the sample by zero airflow.

The analyser uses a 1-hour-cycle analysis with a 40 minute-sampling.

1) Sampling method:




- i. Zero air flow (ml/min)
- ii. Standard gas flow (ml/min)
- iii. Dilution factor: $(i+ii) / ii$

2) Chromatograms:

a) Zero air:

The zero air used for the dilution of the standard bottle must be cleaned at ppt level not to affect the measurement.

- i. Zero air flow : 100ml/min
- ii. Standard gas flow : 0ml/min




Chromatogram of zero air analysis

Our zero air contains some impurities at ppt level:


- Me-cyclohexane: 20-30 ppt
- Ethyl benzene: 20-25 ppt
- N-decane: 10-15 ppt
- 123-TMB: 40-60 ppt
- M-diethylbenzene: 10-15 ppt
- P-diethylbenzene: 10-15 ppt

b) 10ppb:

- i. Zero air flow : 90,2 ml/min
 - ii. Standard gas flow : 10.1 ml/min
 - iii. Dilution factor = 9.9
 - iv. Estimated error on dilution: Relative $\pm 1\%$
- To measure the flows, we use a **Bios, Definer 220-L** flow meter certificated at $\pm 1\%$.




Chromatogram of 10 ppb mixture (area)




Chromatogram of 10 ppb mixture (concentration in ppb)

c) 5ppb:

- i. Zero air flow : 95.0 ml/min
- ii. Standard gas flow : 5.1ml/min
- iii. Dilution factor = 19,6
- iv. Estimated error on dilution: Relative $\pm 1\%$
To measure the flows, we use a **Bios, Definer 220-L** flow meter certificated at $\pm 1\%$.



Chromatogram of 5 ppb mixture (area)




Chromatogram of 5 ppb mixture (concentration in ppb)


d) 1ppb:

- i. Zero air flow : 100,1 ml/min
- ii. Standard gas flow : 0.9 ml/min
- iii. Dilution factor = 111,0
- iv. Estimated error on dilution: Relative $\pm 5\%$

To measure the air flow, we used a **Bios, Definer 220-L** flow meter certificated at $\pm 1\%$. And we use a graduated pipe with a chronometer to measure the flow of standard gas.



Chromatogram of 1 ppb mixture (area)




Chromatogram of 1 ppb mixture (concentration in ppb)


e) 0.1ppb:

- i. Zero air flow : 100.2 ml/min
- ii. Standard gas flow : 0.08 ml/min
- iii. Dilution factor = 1250
- iv. Estimated error on dilution: Relative $\pm$ 5%

To measure the air flow, we used a **Bios, Definer 220-L** flow meter certificated at $\pm 1\%$. And we use a graduated pipe with a chronometer to measure the flow of standard gas.



Chromatogram of 0.1 ppb mixture (area)



Chromatogram of 0.11 ppb mixture (concentration in ppt)

3) Linearity:

We did 4 series of 10 successive analysis with approximate concentration of 10, 5, 1 and 0.1 ppb for each compound.

Three concentrations obtained by the same standard bottle of 99ppb±1ppb diluted by respective factors 9,9/19,6/111/1225.

a) Table of results:


In this table Concentration column represents the concentration obtained by the dilution of benzene (99ppb/dilution factor, see appendix).




Date	Concentration (ppb)	Tube	Volume	2-4-DIME-PENTANE	BENZENE	CYCLOHEXANE	2-3-DIMEC5+2	3-ME-HEXANE	224-TME-PENTANE	N-HEPTANE	ME-CYCLOHEXANE	234-TME-PENTANE	TOLUENE	2-ME-HEPTANE
29/01/2008 23:40	10.0	1	519.85	125582.3	97025.4	90646	211542.2	100401.6	128188	109079.1	110531.6	120949.2	117275	119046.3
30/01/2008 00:40	10.0	2	520.03	125540.6	96660	90534	211334.1	100334.1	128103.1	109039.5	110669.7	120807.2	116501.4	119048.6
30/01/2008 01:40	10.0	3	520.01	125432.1	96256.9	90712	211634.4	100331.6	128096.3	108795.4	110740.6	116709.1		119043.4
30/01/2008 02:40	10.0	4	520.33	125690.2	95969.3	90701.8	211636.4	100452.2	127970.1	109085.8	110891	121087.2	115736.1	119150.8
30/01/2008 03:40	10.0	5	520.51	125800.1	95466.1	90548.4	211308.8	100259.3	127935.6	108702.2	110544.8	120810.1	115030.2	119300.5
30/01/2008 04:40	10.0	6	520.65	125660.5	95286.9	90673.7	211479.8	100288.8	128158.6	109388.3	110526.1	120860.6	114698.1	119199.1
30/01/2008 05:40	10.0	7	520.64	126025.4	95225.6	90729.4	211725.5	100372.8	128369.6	109142.1	110669.1	120981.8	114672.5	119209
30/01/2008 06:40	10.0	8	520.8	125608.2	95139.7	90650.1	211505.1	100271.1	128212.2	109155.8	110294.2	120929	114918.5	119052.2
30/01/2008 07:40	10.0	9	520.97	125752.1	95032.1	90633	211456.7	100130.2	128201.4	108839.6	110556.6	120758	114721.1	119052.8
30/01/2008 08:40	10.0	10	521.06	125651.6	95152.3	90670.7	211210.8	100179.4	128295.4	109151.3	110611.8	120832.2	114666.1	118953.9
30/01/2008 23:02	5.1	1	522.01	62336.9	49495	45033.2	105253.6	50039.9	62036.7	53583	54765.3	59724.5	57626.1	59255.7
31/01/2008 00:02	5.1	2	522.11	62405.5	49106.5	45075.4	105236.5	50020.6	61931.9	53758.3	54692	59781.2	57211.7	59168.9
31/01/2008 01:02	5.1	3	522.04	62859.3	48841.5	45190.2	105694.7	50250.4	62163.6	53867.2	54945.8	59972	56997.3	59435.4
31/01/2008 02:02	5.1	4	521.88	62505.8	48743.6	45276.5	105920.6	50364.2	62240.2	53904.9	55044.2	60127.4	56869.2	59576.9
31/01/2008 03:02	5.1	5	521.88	62864.4	48577.8	45360	106103.2	50333.4	62390.8	54181.6	55249.4	60286.9	56847.8	59724.7
31/01/2008 04:02	5.1	6	521.69	62941	48683.3	45393.9	106174.8	50433.2	62477.9	54612.2	55400.3	60371.7	56844.3	59741.4
31/01/2008 05:02	5.1	7	521.56	63124.5	48805.8	45470	106278	50453.8	62353.3	54123	55351.5	60368.8	56813.2	59732.9
31/01/2008 06:02	5.1	8	521.4	63333.1	48873.2	45595.1	106502.1	50587.5	62702.6	53980.1	55567.6	60700.1	58654.7	59909.8
31/01/2008 07:02	5.1	9	521.36	63263.1	49035.3	45699.2	106713	50714.8	62879.7	54467.2	55339.3	60593.2	57434.1	59928.5
31/01/2008 08:02	5.1	10	521.23	63757.8	49614.7	45904.2	107186.4	50899.7	63220.4	54818.4	55972.6	60935.6	58198.5	60281.2
01/02/2008 18:02	0.9	1	510.04	10579.5	8098.6	7772.5	18535.3	8785.7	10604.2	9079.9	9480.3	9112.1	7480.5	10295.4
01/02/2008 17:02	0.9	2	510.44	10545.6	8090.9	7787.5	18518.2	8780.1	10599.1	9090.5	9398.7	9091.7	7504.8	10236.7
01/02/2008 16:02	0.9	3	510.91	10558.7	8097.4	7720.4	18355.8	8686.8	10625.9	9061.8	9390.6	9022.6	7455.5	10158.4
01/02/2008 19:02	0.9	4	510.27	10535.5	8115	7777.9	18517.7	8742	10721.9	9135	9519.6	9082.1	7437.6	10310.1
01/02/2008 20:02	0.9	5	511.12	10506.8	7984	7693.3	18324.9	8659.1	10563.4	9051.9	9445.9	8981.1	7464.7	10096
01/02/2008 21:02	0.9	6	511.87	10477.2	8015.5	7681.7	18339.3	8551.8	10675.4	8997.4	9452.9	8962.6	7397.5	10129.3
01/02/2008 22:02	0.9	7	512.34	10571.5	7978.3	7698.4	18409.9	8706.4	10655.3	9036.5	9557.5	9026.5	7435.7	10308.3
01/02/2008 23:02	0.9	8	512.7	10572.6	7925.2	7673.3	18370.1	8652.1	10612.4	9073.7	9591.1	9007.8	7427.7	10225.2
02/02/2008 00:02	0.9	9	513.28	10577.7	7945.7	7706.3	18426.7	8667.2	10692.5	9045.8	9694	9034	7437.1	10240.1
02/02/2008 01:02	0.9	10	513.63	10529	8026.6	7776.3	18509.7	8685.1	10641.3	9058.4	9676.1	9071.7	7458.3	10261.8
10/02/2008 03:42	0.08	1	523.1	463.5	587.6	588.9	1552.6	694.4	894.5	729.8	1053	862.2	710.8	832.4
10/02/2008 04:42	0.08	2	523.14	431.5	617.4	619.7	1602.3	714.8	895.8	730.6	1084.5	853.4	694.1	838.6
10/02/2008 05:42	0.08	3	523.24	459.9	661.7	621.7	1541.8	705.2	823.2	697.2	998.4	805.6	686.1	879.4
10/02/2008 06:42	0.08	4	523.33	456.9	625.4	594.1	1568.3	683.3	861.8	744.1	1089.3	833.5	716.9	875.4
10/02/2008 07:42	0.08	5	523.39	473.9	591.2	604.6	1464.7	718.2	860.7	766	1038.2	849.1	709.3	921.6
10/02/2008 08:42	0.08	6	523.66	500.3	622.2	565.7	1661.5	732.3	870.8	726.5	990.1	885.5	723.1	877.1
10/02/2008 09:42	0.08	7	523.72	495.4	649.6	589.4	1673.8	720.2	809.2	747.4	1038.7	860.4	743.3	864.1
10/02/2008 10:42	0.08	8	523.68	475	589.6	609.8	1618	687.4	869.8	768	987.4	834.7	693.9	887.1
10/02/2008 11:42	0.08	9	523.52	451.2	584.7	614.4	1718.3	689.3	779.5	705.5	989.8	847.4	710.9	875.3
10/02/2008 12:42	0.08	10	523.33	469.2	584.6	568.6	1517.9	690.1	800.8	690.9	980.9	859.4	694.1	855.3




Date	Concentration (ppb)	Tube	Volume	3-ME-HEPTANE	N-OCTANE	ETHYLBENZENE	M&P-XYLENES	STYRENE	O-XYLENE	N-NONANE	I-PROPYLENE	N-PROPYLENE	M-ETHYLBENZENE
29/01/2008 23:40	10.0	1	519.85	122378.4	126015.4	105757	233338.3	99322.6	126653.1	131589	122925.3	118005.9	162071.2
30/01/2008 00:40	10.0	2	520.03	122374.8	125119.1	105670.7	233203.9	99280.2	126758.4	131641.7	123094.5	118184.8	163293.4
30/01/2008 01:40	10.0	3	520.01	122474.4	125620.9	105793.3	233352.2	99393.1	126811.1	131660	123140.9	118240.2	163123.8
30/01/2008 02:40	10.0	4	520.33	122474.9	125725.5	105868.5	233537.7	99554.5	127055.5	131808.5	123282.3	118478.1	163303.9
30/01/2008 03:40	10.0	5	520.51	122364.5	125457.6	105608.3	233309	99438.8	126837.1	131736.9	123144	118484.5	163462.2
30/01/2008 04:40	10.0	6	520.65	122501.9	125467.8	105806.7	233375	99619.4	126830.7	131868	123314.4	118602.7	163510.8
30/01/2008 05:40	10.0	7	520.64	122496.4	126083.3	105806.3	233519.6	99671.6	126947.6	131984.8	123392.7	118791.3	162452.7
30/01/2008 06:40	10.0	8	520.8	122673.4	125659.3	105762.9	233256.8	99602.8	126838.5	131816.6	123289.8	118659.2	162161.8
30/01/2008 07:40	10.0	9	520.97	122326.9	125524.8	105571.8	233063.6	99607.6	126838.5	131788.4	123303.5	118574.4	163296.5
30/01/2008 08:40	10.0	10	521.06	122275.1	125504.8	105553.8	232990.2	99664.6	126719.7	131672.9	123202.9	118552.1	163857.1
30/01/2008 23:02	5.1	1	522.01	60916.3	62254.1	53281.9	118289.9	49659.7	61529.4	65635	61545.4	57464.8	79954.2
31/01/2008 00:02	5.1	2	522.11	60668.9	62181.9	53209.5	118037	49520.6	61412.7	65605.6	61535.1	57488	80077.6
31/01/2008 01:02	5.1	3	522.04	61020.3	62328.1	53316.1	118230.2	49653.6	61583.1	65772.8	61743.3	57627.5	80180
31/01/2008 02:02	5.1	4	521.88	61183.2	62521.2	53342.9	118189.4	49740.7	61572.4	65901.4	61831.5	57687.4	80442.9
31/01/2008 03:02	5.1	5	521.88	61259	62550.6	53484.4	118162	49838.5	61749.7	66085.2	61969.7	57837.1	80738.2
31/01/2008 04:02	5.1	6	521.69	61403.7	62619.5	53494.8	118112.1	49863.1	61673.6	66131.1	61969.1	57783.3	80326.4
31/01/2008 05:02	5.1	7	521.56	61354.5	62464.1	53468.7	118099.2	49895.7	61854.9	66203.2	62057.2	57885.5	80958.8
31/01/2008 06:02	5.1	8	521.4	61534.4	62814.9	53591.4	118329	49952	61981	66346.6	62186.1	57961.1	81118.2
31/01/2008 07:02	5.1	9	521.36	61675.3	62923.4	53673.2	118514.1	49988.1	62028.6	66384.2	62246.5	58034.8	80972.4
31/01/2008 08:02	5.1	10	521.23	61989.7	63281.1	53825.5	118936.6	50147.9	62232.3	66590.2	62403.9	58158.5	80837.4
01/02/2008 18:02	0.9	1	510.04	10142.6	10499.6	9452.3	19783.4	7881.4	9171.6	11380.4	10964.6	9030.2	12003.7
01/02/2008 17:02	0.9	2	510.44	10213	10457.5	9383.8	19839.3	7905.2	9135.5	11352.8	11016.1	9020.4	11978.7
01/02/2008 16:02	0.9	3	510.91	10089.7	10373.2	9342.9	19740.1	7863	8858.6	11311.4	10972	9007.2	11982.5
01/02/2008 19:02	0.9	4	510.27	10202.7	10486.7	9344.8	19815	7871.1	9103.1	11276.1	10968.7	9013.4	12018.3
01/02/2008 20:02	0.9	5	511.12	10279.6	10386.1	9393.6	19726.7	7807	9060.3	11387.6	10899.7	8982.3	11929.4
01/02/2008 21:02	0.9	6	511.87	10009.7	10397.7	9309.8	19529.8	7793.9	9013.6	11266.5	10933.3	8931.1	11872.3
01/02/2008 22:02	0.9	7	512.34	10137.7	10460.6	9402	19605.8	7816	9046.2	11311.3	10716.9	8949.1	11931.8
01/02/2008 23:02	0.9	8	512.7	10094.6	10399.5	9521.6	19620.8	7802.1	9058.6	11216.4	10930.4	8938.6	11874.8
02/02/2008 00:02	0.9	9	513.28	10206.2	10385.6	9436.8	19626.5	7810.8	9035.8	11381.1	10846.1	8945.6	11980.9
02/02/2008 01:02	0.9	10	513.63	10286.8	10440.3	9573	19652.4	7814	9091.1	11273.3	10968.1	8978.2	11916.3
10/02/2008 03:42	0.08	1	523.1	882.5	969	1070.8	1774.9	811.7	809.3	1001.4	953.1	901.2	185.1
10/02/2008 04:42	0.08	2	523.14	880.4	938.9	1069.1	1786.4	757.9	851.4	1008.2	939	909.8	192
10/02/2008 05:42	0.08	3	523.24	946	926.2	1056	1775	781.1	910.2	990.8	947.4	935.6	199.2
10/02/2008 06:42	0.08	4	523.33	873.3	935.3	1062.8	1807.4	778.3	855.4	1014	952.2	916.3	179.3
10/02/2008 07:42	0.08	5	523.39	889.5	936	1074.9	1797.9	773.1	878.5	1009.9	946.2	912.8	0
10/02/2008 08:42	0.08	6	523.66	902.6	943.4	1081.9	1794.6	769.8	875.5	1012.9	989.9	909.1	0
10/02/2008 09:42	0.08	7	523.72	938	937.4	1088.5	1786.1	784.7	864.9	1027.7	973.3	911.7	195.6
10/02/2008 10:42	0.08	8	523.68	896.1	937.1	1030.7	1799.4	782.2	863.5	1010	981.7	884.1	191.4
10/02/2008 11:42	0.08	9	523.52	878.9	928	1041.8	1786.5	777.4	859.2	1006.9	949.4	899.6	164.6
10/02/2008 12:42	0.08	10	523.33	868.9	888.3	1066.8	1766.2	771.6	852.5	993.9	962.2	898.1	171

Date	Concentration (ppb)	Tube	Volume	P-ETHYLTOLUENE	135-TMB	O-ETHYLTOLUENE	124-TMB	N-DECANE	123-TMB	M-DIETHYLBENZENE	P-DIETHYLBENZENE	N-UNDECANE
29/01/2008 23:40	10.0	1	519.85	124737.1	124563.5	120929.3	113088	131744.7	104174.1	114325.7	126138.1	117921.1
30/01/2008 00:40	10.0	2	520.03	123563.6	124684.6	120959	113188.5	132022.5	104341.1	114494.9	126313.7	118136
30/01/2008 01:40	10.0	3	520.01	124321.4	124797.6	121195.4	113460.8	132295.8	104532.9	114705.8	126557.9	118204.6
30/01/2008 02:40	10.0	4	520.33	124379.1	125004.5	121421.9	113629.2	132422.7	104665.6	114772.4	126724.8	118257.9
30/01/2008 03:40	10.0	5	520.51	124329.9	124884.9	121327.8	113558.1	132232.6	104650.2	114945.6	126852.1	118444.9
30/01/2008 04:40	10.0	6	520.65	124494	125077	121401.2	113672	132437.4	104792.5	115056.4	126994.3	118681.1
30/01/2008 05:40	10.0	7	520.64	125574.9	125283.9	121559.3	113788	132613.7	104929.8	115212	127239.4	118851.5
30/01/2008 06:40	10.0	8	520.8	125600.5	125192.9	121477.8	113736.4	132493.2	104745.3	115156.3	127252.7	118853.8
30/01/2008 07:40	10.0	9	520.97	124722.4	125110.4	121500.1	113794.4	132522.6	104899.1	115160.1	127340.1	118842.8
30/01/2008 08:40	10.0	10	521.06	123851.8	125168.2	121329.6	113725.4	132409.4	104856.7	115211.4	127273.9	118805
30/01/2008 23:02	5.1	1	522.01	62022.1	62247.4	60901.7	57256.6	66269.9	52987.6	57948.9	63334.8	59795
31/01/2008 00:02	5.1	2	522.11	61890.3	62218.7	60805	57178.8	66211.7	52950.7	57874	63283.7	59696.9
31/01/2008 01:02	5.1	3	522.04	62172.4	62375.4	61009.7	57249	66353.5	53031.7	57925.9	63389.1	59751.6
31/01/2008 02:02	5.1	4	521.88	62045.1	62403.8	61115.5	57348.4	66406.6	53049.7	57987.5	63464.1	59730.5
31/01/2008 03:02	5.1	5	521.88	62040.8	62546.6	61238	57395.1	66568	53171	58143.5	64090.6	59784.4
31/01/2008 04:02	5.1	6	521.69	62615	62549.7	61243.5	57403.5	66575.4	53096.6	58087.4	64064.8	59710.9
31/01/2008 05:02	5.1	7	521.56	61947.1	62573.9	61244.7	57431.6	66601.9	53143.4	58070.9	64065.6	59705.9
31/01/2008 06:02	5.1	8	521.4	62021.8	62666.7	61333.6	57542.5	66741.3	53212.6	58156.9	64205	59796.9
31/01/2008 07:02	5.1	9	521.36	62337.9	62728.8	61421.6	57490	66793.5	53231.5	58269.7	64191.8	59709.7
31/01/2008 08:02	5.1	10	521.23	62967.5	62857.4	61574.8	57683.4	66914.1	53358.8	58324.4	64304.8	59809.1
01/02/2008 18:02	0.9	1	510.04	9391.5	9251.4	10692.8	10122.5	12051.7	9960.7	9712.6	9297.9	10586.2
01/02/2008 17:02	0.9	2	510.44	9422.6	9251.2	10695.4	10141.3	12071.9	9995.3	9735.6	9312	10648.9
01/02/2008 16:02	0.9	3	510.91	9327.1	9204.3	10664.9	10115.7	12114.1	9972.9	9697.1	9265.5	10627.8
01/02/2008 19:02	0.9	4	510.27	9337.4	9261.8	10633.3	10133.6	12017.6	9934.6	9684.3	9285	10591
01/02/2008 20:02	0.9	5	511.12	9335.4	9198.9	10612.7	10106	11992.7	9925	9693.7	9276.4	10597.3
01/02/2008 21:02	0.9	6	511.87	9312.7	9168.5	10520.7	10024.1	11929.1	9860.1	9626.7	9241.2	10693.2
01/02/2008 22:02	0.9	7	512.34	9327.5	9210.9	10521.2	10070.2	11943.2	9860.7	9635.7	9234.5	10523.3
01/02/2008 23:02	0.9	8	512.7	9346.4	9182.1	10523.2	10038.2	11918.4	9859.6	9626	9204.8	10463.1
02/02/2008 00:02	0.9	9	513.28	9235	9199.7	10525	10039.1	11948	9855.4	9605.3	9206.5	10463.6
02/02/2008 01:02	0.9	10	513.63	9364.9	9205.5	10524.3	10051.5	11920.8	9883.3	9622.3	9189.5	10464.8
10/02/2008 03:42	0.08	1	523.1	748.9	897.6	1031.8	928.7	1173.6	1265.8	987.5	939.8	1023.2
10/02/2008 04:42	0.08	2	523.14	751.5	901.2	1054	899.5	1180.4	1306.7	991.8	924.3	1016.1
10/02/2008 05:42	0.08	3	523.24	751.4	915	1028.4	892.3	1174.6	1274	987.9	933.4	1009
10/02/2008 06:42	0.08	4	523.33	747.6	904	1040.7	894.3	1177.4	1259.3	1010.7	937	1013.5
10/02/2008 07:42	0.08	5	523.39	764.4	917.4	1035.9	898.9	1158.3	1265.3	988.9	927.3	1004.5
10/02/2008 08:42	0.08	6	523.66	699.3	904.1	1032.8	885	1170.3	1293.3	1013.2	941.3	972.8
10/02/2008 09:42	0.08	7	523.72	747.8	910.5	1039.6	909.1	1159.4	1272.7	994.1	913.3	998.9
10/02/2008 10:42	0.08	8	523.68	754.4	904	1041.3	906.5	1171.3	1306.4	996.9	928.4	995.4
10/02/2008 11:42	0.08	9	523.52	754.2	913.1	1054.1	889	1124.3	1318.8	1006.8	925.4	1014.1
10/02/2008 12:42	0.08	10	523.33	755.1	890.7	1020.9	896	1182.8	1314.8	1010.1	976.3	1029

b) Linearity Graphs







4) Repeatability:

We calculated the relative standard deviation of areas measured for each concentration and for each compound:

Concentration (ppb)	10,0	5,1	0,9	0,09
	Relative standard deviation (%)			
2-4-DIME-PENTANE	0,13	0,72	0,32	4,34
BENZENE	0,74	0,70	0,86	4,63
CYCLOHEXANE	0,07	0,61	0,58	3,33
2-3-DIMEC5+2MEC6	0,08	0,58	0,45	4,89
3-ME-HEXANE	0,10	0,55	0,79	2,41
224-TME-PENTANE	0,10	0,64	0,45	4,79
N-HEPTANE	0,19	0,73	0,40	3,67
ME-CYCLOHEXANE	0,14	0,70	1,13	4,03
234-TME-PENTANE	0,10	0,65	0,54	2,51
TOLUENE	0,82	1,15	0,40	2,42
2-ME-HEPTANE	0,07	0,56	0,74	2,92
3-ME-HEPTANE	0,09	0,63	0,86	2,96
N-OCTANE	0,22	0,54	0,44	2,13
ETHYLBENZENE	0,10	0,36	0,87	1,65
M&P-XYLENES	0,08	0,22	0,52	0,71
STYRENE	0,14	0,37	0,50	0,78
O-XYLENE	0,09	0,42	0,93	2,95
N-NONANE	0,09	0,50	0,52	1,04
I-PROPYLBENZENE	0,11	0,47	0,78	1,76
N-PROPYLBENZENE	0,20	0,40	0,41	1,50
M-ETHYLTOLUENE	0,37	0,52	0,43	6,62
P-ETHYLTOLUENE	0,52	0,55	0,53	2,36
135-TMB	0,19	0,33	0,34	0,92
O-ETHYLTOLUENE	0,18	0,39	0,72	1,01
124-TMB	0,22	0,26	0,44	1,39
N-DECANE	0,20	0,35	0,58	1,46
123-TMB	0,24	0,24	0,54	1,76
M-DIETHYLBENZENE	0,28	0,26	0,47	1,03
P-DIETHYLBENZENE	0,34	0,65	0,46	1,80
N-UNDECANE	0,30	0,07	0,79	1,59

The relative standard deviation for each series and for each compound is less than 1,2% except for 0.1ppb serie.

5) Conclusion

AirmoVOC C6-C12 has been tested on repeatability (10 successive analysis) and linearity (0.1, 1, 5, 10ppb).

- The determination factors R^2 of linearity graphs are over 0,998.
- For each compound and each concentration, the relative standard deviation of areas is better than 1.2% except for 0.1ppb (6.62% max).

The quality of zero air is very important to have good results with a big dilution (per 1000). An analysis of our zero air (see 2.a) shows some impurities. These impurities increase significantly the concentration measured by the analyser if the sample is a very low concentration (less than 1 ppb). An other problem comes from the integration of certain compounds.

APPENDIX:

Standard Bottle Takachiho certified at $\pm 2\%$

PAMS 58 VOCs CERTIFICATE

N°	Bouteille Takachiho (58 COV)	Formule	PM	BP (°C)	Concentration ppb	Précision (10%) PPB	Précision (+10%) PPB	Concentration µg/m3
1	Ethane	C2H6	30,07	-88	100,8	90,72	110,88	126,08
2	Ethene	C2H4	28,05	-104	102,3	92,07	112,53	119,36
3	Acetylene	C2H2	26,04	-84	101,6	91,44	111,76	110,05
4	Propane	C3H8	44,1	-42,1	100,6	90,54	110,66	184,54
5	Propene	C3H6	42,08	-47,7	104,2	93,78	114,62	182,39
6	Isobutane	i-C4H10	58,12	-12	100,6	90,54	110,66	243,21
7	n-butane	n-C4H10	58,12	-0,5	100,6	90,54	110,66	243,21
8	t-butene	t-2-C4H8	56,11	1	102,0	91,8	112,2	238,07
9	1-butene	1-C4H8	56,11	-6,3	103,2	92,88	113,52	240,87
10	c-2-butene	c-2-C4H8	56,11	3,7	100,1	90,09	110,11	233,64
11	isopentane	i-C5H12	72,15	30	100,1	90,09	110,11	300,42
12	1-pentene	1-C5H10	70,14	29,9-30,1	100,6	90,54	110,66	293,51
13	n-pentane	n-C5H12	72,15	35-36	100,6	90,54	110,66	301,93
14	t-2-pentene	t-2-C5H10	70,14	37	118,4	106,56	130,24	345,45
15	cyclopentane	(CH2)5	70,14	50	98,5	88,65	108,35	287,39
16	c-2-pentene	c-2-C5H10	70,14	37-38	94,4	84,96	103,84	275,42
17	2-methyl-1,3-butadiene	2-(CH3)-1,3-C4H5	68,12	34	100,3	90,27	110,33	284,21
18	2,2-dimethylbutane	2,2-(CH3)2C4H8	84,16	50	98,3	88,47	108,13	344,13
19	2,3-dimethylbutane	2,3-(CH3)2C4H8	86,18	58	97,3	87,57	107,03	348,81
20	2-methylpentane	2-(CH3)C5H11	86,18	62	97,6	87,84	107,36	349,88
21	3-methylpentane	3-(CH3)C5H11	86,18	64	97,5	87,75	107,25	349,52
22	2-methyl-1-pentene	2-(CH3)-1-C5H9	84,16	62	97,1	87,39	106,81	339,93
23	n-hexane	n-C6H14	86,18	69	98,9	89,01	108,79	354,54
24	methylcyclopentane	(CH3)C5H9	84,16	72	97,3	87,57	107,03	340,63
25	2,4-dimethylpentane	2,4-(CH3)2C5H10	100,21	80	98,3	88,47	108,13	409,76
26	benzene	C6H6	78,11	80	99,0	89,1	108,9	321,67
27	cyclohexane	(CH2)6	84,16	80,7	97,6	87,84	107,36	341,68
28	2-methylhexane	2-(CH3)C6H13	100,21	90	97,5	87,75	107,25	406,43
29	2,3-dimethylpentane	2,3-(CH3)2C5H10	100,21	89-90	98,4	88,56	108,24	410,18
30	3-methylhexane	3-(CH3)C6H13	100,21	91	97,6	87,84	107,36	406,84
31	2,2,4-trimethylpentane	2,2,4-(CH3)3C5H9	114,23	98-99	97,5	87,75	107,25	463,29
32	n-heptane	n-C7H16	100,21	98	97,7	87,93	107,47	407,26
33	methylcyclohexane	(CH3)C6H11	98,19	101	100,1	90,09	110,11	408,85
34	2,3,4-trimethylpentane	2,3,4-(CH3)3C5H9	114,23	113-114	98,5	88,65	108,35	468,04
35	toluene	(CH3)C6H5	92,14	110,6	97,4	87,66	107,14	373,31
36	2-methylheptane	2-(CH3)C7H15	114,23	116	98,0	88,2	107,8	465,66
37	3-methylheptane	3-(CH3)C7H15	114,23	118	98,3	88,47	108,13	467,09
38	n-octane	n-C8H18	114,23	125-127	98,6	88,74	108,46	468,51
39	ethylbenzene	(C2H5)C6H5	106,17	136	97,4	87,66	107,14	430,16
40	m-xylene	m-C8H10	106,17	138-139	100,3	90,27	110,33	442,96
41	p-xylene	p-C8H10	106,17	138	97,3	87,57	107,03	429,71
42	styrene	(C2H3)C6H5	104,15	145-146	98,5	88,65	108,35	426,74
43	o-xylene	o-C8H10	106,17	143-145	99,8	89,82	109,78	440,76
44	n-nonane	n-C9H20	128,26	151	98,1	88,29	107,91	523,39
45	isopropylbenzene	i-(C3H7)C6H5	120,2	152-154	100,2	90,18	110,22	501,00
46	α-pinene	α-C10H16	136,24	155-156	99,4	89,46	109,34	563,32
47	n-propylbenzene	n-(C3H7)C6H5	120,2	159	100,0	90	110	500,00
48	m-ethyltoluene (3-)	m-(C2H5)C7H7	120,2	158-159	100,2	90,18	110,22	501,00
49	p-ethyltoluene (4-)	p-(C2H5)C7H7	120,2	162	99,7	89,73	109,67	498,50
50	1,3,5-trimethylbenzene	1,3,5-(CH3)3C6H3	120,2	162-164	100,3	90,27	110,33	501,50
51	o-ethyltoluene (2-)	o-(C2H5)C7H7	120,2	164-165	99,9	89,91	109,89	499,50
52	β-pinene	β-C10H16	136,24	164-165	99,5	89,55	109,45	563,89
53	1,2,4-trimethylbenzene	1,2,4-(CH3)3C6H3	120,2	168	100,2	90,18	110,22	501,00
54	decane	n-C10H22	142,29	174	99,1	89,19	109,01	586,56
55	1,2,3-trimethylbenzene	1,2,3-(CH3)3C6H3	120,2	175-176	100,1	90,09	110,11	500,50
56	m-diethylbenzene (1,3-)	m-(C2H5)2C6H4	134,22	183	100,2	90,18	110,22	559,44
57	p-diethylbenzene (1,4-)	p-(C2H5)2C6H4	134,22	184	99,3	89,37	109,23	554,41
58	undecane	n-C11H24	156,31	196	100,0	90	110	650,21
59	Azote	N2			BALANCE			

At 20°C under 1atm: Molar Weight = 24.04 L /mol