


DRIVER MODBUS/JBUS TRAINING

Chromatotec®

Overview



COMMUNICATION PROTOCOLS

- Considering PC networks and telecommunications, a ***communication protocol*** is defined to be a specification of several rules dedicated to a specific communication type.
 - Initially, a protocol was defined as a set of rules used to communicate on a same « abstraction level » between two different machines. By extension we also use today this word to talk about the communication rules between two abstraction levels on a same instrument.
- The mostly used protocol is the network one.

DEFINITION : MODBUS

- **Modbus** and **Jbus** are *communication protocols* used for *programmable automations* networks. They are working on Master / Slave mode. They are made up of frames containing the selected automaton address, the action to be done, function, (writing, reading), the data and the error check code called CCR (control of cyclic redundancy).
- The *Modbus* protocol (**Modicon trade mark**) and the *Jbus* one (**April trade mark**) are dialogue protocols based on a hierarchical structure between a master and several slaves.



OVERVIEW

There are two ways to communicate :

1. The Master talks to the Slave and waits for an answer.
2. The master talks to all the slaves without waiting any answer (general broadcasting)


The Master is able to address up to 255 slaves.

COMMUNICATION RULES


- The Master handles alone the exchanges : **He decides**,
- Repeat the question if an erroneous exchange,
- Consider a slave missing after a time-out,
- A single equipment at the same time in emission, on one line,
- No slave can send message without request from the master,
- Possible side communication between slaves if and only if the Master was programmed to receive data and send them from a slave to an otherone.

TREATED ON A HIERARCHICAL BASIS STRUCTURE

One Master, several slaves



FRAME JBUS / MODBUS



Slave # : 0 to 255 (0 = general broadcast for writing functions only).

Function Code : Indicates the exchange type (see next pages)


Data Field : Information field containing the parameters linked to the function : word address, word value, number of words...

CRC16 : Word used to detect transmission errors (Cyclic Redundary Check).

FRAME : The Query ?

It contains a Function Code indicating to the sent addressed slave, which type of action is asked. The data contain additional information which the slave needs to execute this function.

The CRC16 allows the slave to make sure of the completeness of the contents of the question.



FRAME : The Response

Slave Address 11


Function 04

Byte Count 02

Data Hi (Register 30009) 00

Data Lo (Register 30009) 0A

Error Check (LRC or CRC) —




FRAME : The Response

If an error appears, function code is modified to indicate that the answer is an error answer.

The data contain then an exception code identifying the error type.

The control field allows the master to confirm that the message is valid.



FRAME : The coding

Two types of code can be used to communicate on a Modbus network. All the present equipments on the network must be configured according to the same type.

ASCII Type : each component of the frame is coded with ASCII characters (2 x 8 bits)

RTU Type (Remote Terminal Unit) : each byte composing a frame is coded with hexadecimal characters (2 x 4 bits). The maximal size of the data is 256 bytes.

The RTU mode allows a higher data flow for the same transmission speed.

MODBUS FUNCTIONS

The MODBUS protocol is composed of predefined functions, associated with a particular code, which we can be classified in 3 families :

1.- Data access : these functions allow to reach in reading and/or writing mode, bits, words or files of ModBus equipments.

2.- Diagnostics : these functions allow to make some diagnose on a ModBus equipment.

3.- Others : These functions allow the encapsulation of the MODBUS protocol on another one (ex: Can open)




MJBus CHROMATOSUD Driver Overview

- ✓ The MJBus driver is a program that realizes the interface between the software Vistachrom, that pilots the chromatograph and the ModBus (RTU mode) or Jbus communication protocol. This driver uses the Vistachrom real time database (RTDB or BDTR : Base de Donnée Temps Réel, in french).

- ✓ This is a shared memory area that contains amongst other things, the latest results received from the instruments.

MJBus CHROMATOSUD Driver Overview





MJBus CHROMATOSUD Driver Overview

- ✓ This driver can run several slave numbers if necessary. Currently it implements, the following 3 functions of the protocol:
 - Reading of n words.
 - Writing of a word
 - Simultaneous writing of several words


MJBus CHROMATOSUD Driver Configuration : Mapping

✓ Before using the MJBus driver, the configuration file must be checked and edited. To do it you have to:

- You should have a configured file, let's see what it looks like. The configuration file will define what we call the MAPPING of MJBus driver.




MJBus CHROMATOSUD Driver Overview



MJBus CHROMATOSUD Driver Configuration : Defines

✓ These details are divided in several sections. These sections begin with a name of section, this name is framed with brackets.

• The first section called [Defines] allows to set few functions and to define the symbolic constant that will be used in the other sections in order to make the writing or the configuration easier.



```
DRVJMJBUS_145_Model_ModBus.INI - Bloc-notes
Fichier Edition Format Affichage ?
*****
*
*                               SETUP MODBUS
*                               Analyzer #nnnnnnn
*       March 2007
*****
[Defines]

RateCycleTimeout=150
HoldError=60

$Instrument1=#nnnnnnn
$Sequence1=""

$Measure=(Normal)
$Calibration=(Calib)
$Zero=(Zero)

$Defaultvalue=65535
$ExpScaleConcentration=2
$ExpScaleRetentionTime=1
$ExpScaleAreaPic=-1

$Substance1="$Instrument1,$Sequence1,$Measure,(Subst1)"
$Substance2="$Instrument1,$Sequence1,$Measure,(Subst2)"
$Substance3="$Instrument1,$Sequence1,$Measure,(Subst3)"
$Substance4="$Instrument1,$Sequence1,$Measure,(Subst4)"
$Substance5="$Instrument1,$Sequence1,$Measure,(Subst5)"
$Substance6="$Instrument1,$Sequence1,$Measure,(Subst6)"
$Substance7="$Instrument1,$Sequence1,$Measure,(Subst7)"
$Substance8="$Instrument1,$Sequence1,$Measure,(Subst8)"
$Substance9="$Instrument1,$Sequence1,$Measure,(Subst9)"

$Substance1Calib="$Instrument1,$Sequence1,$Calibration,(Subst1)"
$Substance2Calib="$Instrument1,$Sequence1,$Calibration,(Subst2)"
$Substance3Calib="$Instrument1,$Sequence1,$Calibration,(Subst3)"
$Substance4Calib="$Instrument1,$Sequence1,$Calibration,(Subst4)"
$Substance5Calib="$Instrument1,$Sequence1,$Calibration,(Subst5)"
```



MJBus CHROMATOSUD Driver Configuration : Defines


✓ SECTION CONFIGURATION [DEFINES] : This section can modify the value of two parameters and to define the characters string symbolic constants.

- RATE CYCLE TIME OUT : This parameter indicates the maximum time accepted between two measures series in rate cycle time out (percentage). It defines the maximum allowed delay for the analyser to release its next results.
- DEFINITION OF SYMBOLIC CONSTANTS : In order to make the writing and configuration of definition sections of the Modbus slave mapping ([Slave n]) easy, we can define the symbolic constants that stand for a string.
- Save time!

MJBus CHROMATOSUD Driver

Configuration : Slaves overview

- The following sections allow to define the ModBus / Jbus registers « mapping » of the different slaves necessary to the setting up. The name of these sections is in the form [Slave n] where n is the ModBus/Jbus slave number. You can define as many slave as necessary for you configuration but generally one slave is sufficient.

A screenshot of a Windows-style application window titled "Edit MJBus Driver 1 setup file: : D:\WISTAC~1\MJBUS_~1\DRV\MJBUS_144_1.INI". The window contains a text editor with the following configuration for "Slave 1":

```
File
;===== ModBUs Mapping Slave 1 =====
[Slave 1]
0=wSamplingDay,$Instrument1,$Sequencel,$Measure
1=wSamplingMonth,$Instrument1,$Sequencel,$Measure
2=wSamplingYear,$Instrument1,$Sequencel,$Measure
3=wSamplingHour,$Instrument1,$Sequencel,$Measure
4=wSamplingMinute,$Instrument1,$Sequencel,$Measure
5=wSamplingSecond,$Instrument1,$Sequencel,$Measure

6=wSamplingDay,$Instrument1,$Sequencel,$Zero
7=wSamplingMonth,$Instrument1,$Sequencel,$Zero
8=wSamplingYear,$Instrument1,$Sequencel,$Zero
9=wSamplingHour,$Instrument1,$Sequencel,$Zero
10=wSamplingMinute,$Instrument1,$Sequencel,$Zero
11=wSamplingSecond,$Instrument1,$Sequencel,$Zero

12=wSamplingDay,$Instrument1,$Sequencel,$Calibration
13=wSamplingMonth,$Instrument1,$Sequencel,$Calibration
14=wSamplingYear,$Instrument1,$Sequencel,$Calibration
```

At the bottom of the window, there are two buttons: "Save & Quit" and "Quit".

MJBus CHROMATOSUD Driver

Configuration : Slaves Mapping definition1

- ***Important note*** : the register numbering varies between Modbus and Jbus protocol, despite the other parts are identical. The Modbus protocol starts the register numbering with 1 (or 40001) and the Jbus starts with 0. Depending on the Modbus client, this may lead that the register address may be 1 shifted and one should take in account when building request.
- Each line of this section defines a field in the slave register space. According to its value type, a field can use one or more registers to fit.

MJBus CHROMATOSUD Driver Configuration : Slaves Mapping definition2

- The syntax of these lines is the following : Address=FieldType, parm1, parm2...

- Example :

```
100=wConcentration, #5411002, Seq541, Mth541_1, Toluene,  
2, 65535
```

- Assigns to register address 100, the integer of the toluene. And this on the condition that the concentration has been obtained with the Mth541_1 method of the Seq541 sequence running on the #5411002 instrument. The default value, without any measurement in the allotted time, is 65535.

MJBus CHROMATOSUD Driver

Configuration : Slaves Mapping definition3

- LIST OF DIFFERENT FIELD TYPES

| File type | Description | Size | Type |
|------------------------|---|-------------|-------------|
| sConstant | Inserts a fix text | - | String |
| sDate | Inserts the PC date (in the form of a string) | 4 | String |
| sTime | Inserts the PC time (in the form of a string) | 3 | String |
| sSamplingDate | Inserts the date of the latest sampling achieved with a data method (in the form of a string) . | 4 | String |
| sSamplingTime | Inserts the hour of the latest sampling achieved with a data method (in the form of a string) . | 3 | String |
| wSamplingDay | Inserts the date of the latest sampling achieved with a data (in the form of a integer) . | 1 | Word |
| wSamplingMonth | Inserts the month of the latest sampling achieved with a data method (in the form of a integer). | 1 | Word |
| wSamplingYear | Inserts the year of the latest sampling achieved with a data method (in the form of a integer) . | 1 | Word |
| wSamplingHour | Inserts the hour of the latest sampling achieved with a data method (in the form of a integer). | 1 | Word |
| wSamplingMinute | Inserts the minutes of the latest sampling achieved with a data method (in the form of a integer).. | 1 | Word |
| wSamplingSecond | Inserts the seconds of the latest sampling achieved with a data method (in the form of a integer). | 1 | Word |
| fSamplingDay | Inserts the day of the latest sampling achieved with a data method (in the form of a float) . | 2 | Float |
| fSamplingMonth | Inserts the month of the latest sampling achieved with a data method (in the form of a float) . | 2 | Float |
| fSamplingYear | Inserts the year of the latest sampling achieved with a data method (in the form of a float) | 2 | Float |

MJBus CHROMATOSUD Driver

Configuration : Slaves Mapping definition4

- LIST OF DIFFERENT FIELD TYPES

| | | | |
|------------------------|---|---|---------------------|
| fSamplingYear | Inserts the year of the latest sampling achieved with a data method (in the form of a float) | 2 | Float |
| fSamplingHour | Inserts the hour of the latest sampling achieved with a data method (in the form of a float). | 2 | Float |
| fSamplingMinute | Inserts the minutes of the latest sampling achieved with a data method (in the form of a float). . | 2 | Float |
| fSamplingSecond | Inserts the seconds of the latest sampling achieved with a data method (in the form of a float) . | 2 | Float |
| wConcentration | Inserts the concentration of a substance (in the form of a integer) | 1 | Word |
| wRetentionTime | Inserts the retention time of a substance (in the form of a integer) | 1 | Word |
| wPicArea | Inserts the peak surface of a substance (in the form of a integer). | 1 | Word |
| lConcentration | Inserts the substance concentration (in the form of a integer) | 2 | Long |
| lRetentionTime | Inserts the retention time of a substance (in the form of a long integer) | 2 | Long |
| lPicArea | Inserts the peak surface of a substance (in the form of a long integer) | 2 | Long |
| fConcentration | Inserts the substance concentration (in the form of a float) | 2 | Float |
| fRetentionTime | Inserts the retention time of a substance (in the form of a float) | 2 | Float |
| fPicArea | Inserts the surface of a substance peak (in the form of a float) | 2 | Float |
| a2Concentration | Inserts the substance concentration (in the form of a float number with 2 characters identifying the substance before) | 3 | CodeSubst +Float |
| a2RetentionTime | Inserts the retention time of a substance (in the form of a decimal float number with two characters identifying the substance before) | 3 | CodeSubst +Float |
| a2PicArea | Inserts the surface of the substance peak (in the form of a float number with two characters identifying the substance before). | 3 | CodeSubst +Float |

MJBus CHROMATOSUD Driver

Configuration : Slaves Mapping definition5

- LIST OF DIFFERENT FIELD TYPES

| | | | |
|----------------------|--|---|----------|
| Command | Inserts a order register allowing to act on the instrument (log on/log off the cycle, calibration, zero) | 2 | Reg. 32b |
| Status | Inserts a state register reporting the instrument state (log on/log off, run/standby) | 4 | Reg. 64b |
| Default | Inserts a default register reporting some error codes emitted by the instrument | 2 | Reg. 32b |
| word | Inserts a RTDB data (non signed integer of 16 bits) | 1 | Word |
| int | Inserts a RTDB data (signed integer of 16 bits) | 1 | Int |
| long | Inserts a RTDB data (signed integer of 32 bits) | 2 | Long |
| float | Inserts a RTDB data (floating comma integer) | 2 | Float |
| wLifeSignal | Inserts the « live signal » of an instrument.. In the RTDB it means to the « LifeSignal » rubric value of the instrument key (ex : #5000404.LifeSignal). This value is increased roughly one time by second.. | 1 | Word |
| wResultsCount | Inserts the « results counter » of an instrument. In the.RTDB, it means to the « ResultsCount » rubric value of the instrument key (ex : #5000404.ResultsCount). This value is increased each time when the instrument supplies new results.. | 1 | Word |

MJBus CHROMATOSUD Driver

Declaration of Instrument register fields1

Command register

| Command N° | Name | Action |
|------------|--------|--|
| 1 | LogOn | Asks to Vistachrom to connect to the instrument |
| 2 | LogOff | Asks to Vistachrom to disconnect to the instrument |
| 3 | Start | Asks to Vistachrom to start up the sequence |
| 4 | Stop | Asks to Vistachrom to stop the current sequences |
| 5 | Calib | Asks to Vistachrom to insert a calibration method |
| 6 | Zero | Asks to Vistachrom to insert a zero method |
| 7 | Mth1 | Asks to Vistachrom to insert the method 1 |
| 8 | Mth2 | Asks to Vistachrom to insert the method 2 |
| 9 | Mth3 | Asks to Vistachrom to insert the method 3 |
| 10 | Mth4 | Asks to Vistachrom to insert the method 4 |
| 11 | Mth5 | Asks to Vistachrom to insert the method 5 |
| 12 | Mth6 | Asks to Vistachrom to insert the method 6 |

MJBus CHROMATOSUD Driver

Declaration of Instrument register fields2

State Register

| Bit number of the internal state | Meaning |
|----------------------------------|---|
| 33 | Log On -> The instrument is connected to Vistachrom |
| 34 | Run -> The instrument executes a sequence |
| 35 | Wait -> The instrument is on the waiting phase (above all useful when it is in the « slave » mode) |
| 36 | Synch -> th instrument is in progress of synchronisation (ex : waiting for the whole minute for starting up a new method)) |
| 37 | Calib -> The latest results obtaines come from the calibration method. |
| 38 | Zero -> The latest results obtained come from the zero method. |
| 39 | ResultMth1 -> the latest obtained results come from the method 1 |
| 40 | ResultMth2 -> The latest obtained results come from the method 2 |
| 41 | ResultMth3 -> The latest obtained results come from the method 3 |
| 42 | ResultMth4 -> The latest obtained results come from the method 4 |
| 43 | ResultMth5 -> The latest obtained results come from the method 5 |
| 44 | ResultMth6 -> The latest obtained results come from the method 6 |
| 45 | SamplingMth1 -> the sampling phase of the method 1 is pending |
| 46 | SamplingMth2 -> the sampling phase of the method 2 is pending |
| 47 | SamplingMth3 -> The sampling phase of the method 3 is pending |
| 48 | SamplingMth4 -> The sampling phase of the mehtod 4 is pending |
| 49 | SamplingMth5 -> The sampling phase of the method 5 is pending |
| 50 | SamplingMth6 -> The sampling phase of the method 6 is pending |
| 51 | TimeOut -> the insturment has not returned any results in the given time. |
| 52 | Reserved |
| | |
| 64 | Reserved |


MJBus CHROMATOSUD Driver

INTERFACE DRIVER USER1

THE « SETUP & LOG » TAB

THE “ SETUP COM PORT ” FRAME

It allows to setup the asynchronous serial communication port.



MJBus CHROMATOSUD Driver


INTERFACE DRIVER USER2

THE « SETUP & LOG » TAB

THE FRAME « LOG »

The display area shows the actions executed by the driver.

- The configuration file reading (DrvMJBus_14x_1.INI)
- The initialization of RateCycleTimeOut
- The displayed errors occurred when the configuration file is interpreted.



MJBus CHROMATOSUD Driver

INTERFACE DRIVER USER4


THE « MAPPING » TAB

This tab visualizes the « mapping » of different ModBus / JBus slaves run by the driver. This mapping is under the form of a table. Each line of this one represents a field in the slave register space

MJBus CHROMATOSUD Driver

INTERFACE DRIVER USER4

THE « MAPPING » TAB



Green: indicates that the driver could open the serial port.

Indicates the configuration files to be used.

Indicates the access level

| Address | Label | Value | Type | Size | Source parameter |
|---------|-----------------|------------------|--------|------|----------------------|
| 5 | Manufacturer | Chromatotec | String | 10 | |
| 10 | Reference | BTX1000 | String | 5 | |
| 15 | Identification | #2201006 | String | 5 | |
| 20 | Soft Version | Vistachrom 1.4.4 | String | 10 | |
| 30 | System Date | 22/02/2007 | String | 4 | |
| 34 | System Time | 18:26:31 | String | 3 | |
| 37 | Sampling Day | 30 | Float | 2 | #2201006, % (Normal) |
| 39 | Sampling Month | 12 | Float | 2 | #2201006, % (Normal) |
| 41 | Sampling Year | 1899 | Float | 2 | #2201006, % (Normal) |
| 43 | Sampling Hour | 0 | Float | 2 | #2201006, % (Normal) |
| 45 | Sampling Minute | 0 | Float | 2 | #2201006, % (Normal) |
| 47 | Sampling Second | 0 | Float | 2 | #2201006, % (Normal) |

MJBus CHROMATOSUD Driver

RUNNING MJBUS DRIVER

- Now let's practice : hands on the GCs

