

# Le Bulletin

## CONTENTS

|                                  |    |
|----------------------------------|----|
| Our history                      | p1 |
| Certifications and references    | p1 |
| Certifications of our equipments | p1 |
| Training distributors/customers  | p2 |
| Our team                         | p2 |

## EXHIBITIONS 2010

### POLLUTEC 2009 - FRANCE

Paris – December 1<sup>st</sup> to 4<sup>th</sup> 2009

<http://www.pollutec.com/>

### ARAB LAB 2010 – UAE

Dubai – January 9<sup>th</sup> to 12<sup>th</sup> 2010

<http://www.arablab.com/home.html>

### ANALYSE INDUSTRIELLE

Paris – February 2<sup>nd</sup> to 4<sup>th</sup> 2010

### ACHEMASIA 2010 – CHINE

Beijing – June 1<sup>st</sup> to 4<sup>th</sup> 2010

<http://www.achemasia.de>

### AWMA 2010 – CANADA

Calgary - June 22<sup>nd</sup> to 25<sup>th</sup> 2010

<http://www.awma.org/ACE2010/>

## Our history

Since its creation in 1986 the company Chromato-Sud has continuously expanded.

In 1992, Chromato-Sud created its own gas analyzers.

In 1996, the Chromatotec Group was born by integrating into Chromato-Sud (France) the companies airmotec GmbH (Germany), airmotec AG (Switzerland) RN Electronique (France).

In 2001, we developed our own software for the acquisition and treatment of data : **VISTACHROM**.

Our instruments have kept on evolving in order to meet our customers' requirements.

From 2001, Chromatotec has started supplying complete analytical cabinets for various applications.

Thus, to date, we have sold more than 1000 "turn-key" solutions worldwide.

## Certifications and references

Our analyzers, which are at the cutting edge of technology, have gained specific certifications relating to the relevant standards. Furthermore these tests were performed by bodies recognised worldwide (for example US EPA).

### USA :

ASTM norm for Energy Medor (20 Tests US EPA for airTOXIC in 2009

### Europe :

TUV and CNR certifications (1996 and 2006) Mcerts certification underway with NPL (2010)  
For models : airmoVOC C6C12 and airTOXIC

### China :

Beijing institute and Shanghai institute

Combined with our continuous product development, this makes us 'The Experts in Gas Analysis'.



## Certifications of our instruments

Our Medor analyzers (on-line detection of sulfurs in air or natural gas/ LPG) have received the ISO 6326/2 norm in 1981 : gas chromatographic method for the qualitative and quantitative analysis of odoriferous sulfur compounds in natural gas.

Our energyMEDOR analyzer has been chosen as the reference instrument the standard method ASTM D7493-08 since April 2009.



Our airmoBTX1000 instrument obtained new certification from the Italian CNR in July 2007.



Our airTOXIC analyzer was certified in 2006 by an approved European organization, CNR, for the detection of BTEX by PID detector.

Our airmoBTX1000 was certified in 1996 by the German TÜV organization. The study done on our analyzer was the analysis of benzene, toluene, ethyl benzene, m-p xylene, x-o Xylene in air monitoring according to the national standard method.



Our latest instrument Medor exp dedicated to the measurement of sulfur compounds in natural gas and fuels has been approved by CSA in June 2009. Specification Ex Class I Div 2, group C&D.



## Trainings for distributors / customers

Chromatotec regularly organizes training sessions for its distributors / customers within its facilities of Saint-Antoine (near Bordeaux, France). During these sessions, visitors are trained on the use of our instruments.

Recently, we have organized two training sessions for distributors in early September and early October.

The topics of these trainings sessions were :

- Understanding of our analyzer schematics and parameters
- Easy start set-up and calibration of the analyzers
- Training on Vistachrom software
- Training on preventive maintenance and trouble shooting

Next training sessions for distributors will take place on :

- December 7<sup>th</sup> and 8<sup>th</sup> 2009, we will receive German and Chinese distributors for training on airMOzone and airTOXIC.
- January 13<sup>th</sup> to 15<sup>th</sup> 2010

In addition to training we organize specific one-day seminars on more general themes, such as :

- Seminar on VOC which took place at our facilities in January 2009.
- Seminar on Sulfurs : in march 2009 at our facilities and in October 2009 in Paris.

DVDs for maintenance are now available :

- Service TRSMEDOR
- Service energymEDOR
- Service airTOXIC
- Service airMOzone
- Service chromaS
- Service PID & FID analyzers

Contact our Customer service at: [support@chromatotec.com](mailto:support@chromatotec.com) to receive one of these DVDs providing the serial number of the analyzer.



## Sales team [info@chromatotec.com](mailto:info@chromatotec.com)



Régis FRANCOIS  
Sales manager France  
[regis.francois@chromatotec.com](mailto:regis.francois@chromatotec.com)



Franck AMIET  
President & General manager  
Sales manager China / Korea / Japan / South America  
[franck.amiet@chromatotec.com](mailto:franck.amiet@chromatotec.com)



Mike BOWEN  
General EXPORT Sales manager  
Sales manager  
Australia / Central Asia / Middle East / India / Russia  
[mike.bowen@chromatotec.com](mailto:mike.bowen@chromatotec.com)



Yan PASCOLI  
Executive & Sales assistant  
[yan.pascoli@chromatotec.com](mailto:yan.pascoli@chromatotec.com)

Anne-Solène PALCY  
Sales assistant  
[anne-solene.palcy@chromatotec.com](mailto:anne-solene.palcy@chromatotec.com)

Marie-Pierre MAASS  
Sales assistant  
[marie-pierre.maass@chromatotec.com](mailto:marie-pierre.maass@chromatotec.com)



Michel ROBERT  
Customer service manager  
[michel.robert@chromatotec.com](mailto:michel.robert@chromatotec.com)



Bertrand JULY  
Customer service engineer  
[bertrand.july@chromatotec.com](mailto:bertrand.july@chromatotec.com)



Benjamin LAMARQUE  
Customer service engineer  
[benjamin.lamarque@chromatotec.com](mailto:benjamin.lamarque@chromatotec.com)



Didier WANG XIAOMING  
Marketing and sales support (China)  
[didierwxm@163.com](mailto:didierwxm@163.com)

[support@chromatotec.com](mailto:support@chromatotec.com)  
**Service team**

## CHROMATOTEC

**AMERIQUE – UNITED STATES**  
18333 Egret Bay Blvd, Suite 270  
Houston TX 77058  
Ph +1 281 335 4944  
Fax +1 281 335 4943



**EUROPE – FRANCE**  
15 Rue d'Artiguelongue  
33240 Saint-Antoine  
Ph. +33 (0)557 940 626  
Fax +33 (0)557 940 620



**ASIE – CHINE**  
Room 765, Office Tower, Poly Plaza  
No.14, Dongzhimen Nandajie,  
Dongcheng District  
Beijing 100027  
Ph. + 86 10 65 52 75 98  
Fax + 86 10 65 53 38 98



[info@chromatotec.com](mailto:info@chromatotec.com)