

# airmoTWA分析系统机柜

自动连续监测化学组分  
Automatic Continuous Monitoring of chemical compounds

GC/MS冷阱  
Trap GC / MS

6通道多路进样器  
Multiplexer 6 ways

6通道多路进样器 + 校准仪  
Multiplexer 6 ways + calibrator

监视器  
Supervisor

带其加热传输管道的质谱仪:

DET QMS  
Mass Spectrometer with heated transfer line:  
DET QMS

airmo VOC Expert气相色谱分析仪  
airmoVOC Expert

氢气发生器: hydroxichrom  
H2 generator: Hydroxichrom

监测以下组分:

- \* T014 / T015 / PAMS表组分
- \* 有毒组分 (有机金属化合物、有机硅烷)
- \* 水中VOC组分 (502-2 / 524 方法)

Surveillance of:

- T014 / T015 / PAMS
- Toxic compounds (organometallics, organosilanes)
- Volatiles compounds from water (502-2 / 524)

各种应用:

- 洁净室
- 工艺过程
- 工业气体

Different applications:

- Clean Air Room
- Process
- Industrial air

